

# NIBP Tester / Simulator

# SimCube



The **SimCube** is an NIBP Simulator designed for evaluating and testing NIBP monitors.

The hand-held **SimCube** is the smallest and easiest to use NIBP monitor tester available. Operation is via one push button. Easily select one of four built-in NIBP test modes. Add the ECG, Respiration, and Pacer simulation option to evaluate patient monitors. Test values are displayed on a large LED display.

**SimCube** - NIBP Tester is the instrument of choice, for routinely testing and maintaining NIBP monitors.

## Features:

- **Rugged, Accurate, Best Value**
- **Portable, Easy to Use. One Button Operation.**
- **Patient Simulation - ECG - Respiration - Pacer NIBP Simulation**
- **Modes Manometer, Adult, Neonatal and Hypertensive**

## Specifications:

**Display:** LED 3 ½ digits

**Controls:** Push button Mode selector LED Indication of mode

**Connector:** Female Luer lock connector

**Operating Voltage:** 110 to 240 VAC Adapter

**Dimensions:** 3 x 3 x 3.5 inches( 7.6 x 7.6 x 8.9 cm)

**Weight:**2.5 lbs (1.1. kg) Measurements:

**Manometer:** 0 to 480 mm Hg

**Accuracy:** ± 1% of reading

**Adult:** 120/80 (100) mmHg 69 BPM Heart Rate:

**Neonatal:** 70/40 (55) mmHg 80 BPM

**Hypertensive:** 190/120 (150) mmHg 69 BPM

**ECG:** 1 mV Amp R Wave. Respiration: 4 W SQ Wave 35 BPM, Pacer: 3 mV Amp, 1.2 ms width

**For Product updates and accessories**

**Please visit**

<http://NIBPSimulator.com>

## Ordering Information:

Part #	Description	Part #	Description
VSC-1	NIBP Tester	VADAPT-D	Dynamp Adapter
VSC-2	NIBP Simulator Include VAPAT-D, Q, L or M	VADAPT-Q	Quick Connect Adapter
VSC-CASE	Carrying Case	VADAPT-L	Luer Lock Connect Adapter
		VADAPT-M	Marquette Adapter



60 Bethpage Drive, Hicksville, NY 11810, Phone: 516-433-7400, Fax: 516-433-7458

E-Mail [Sales@Netech.org](mailto:Sales@Netech.org)

<http://www.GoNetech.com>